

Supplementary Table 1. Antibodies and reagents used in this study

Antibody/reagent	Clone	Company	Dilution
CD3	Polyclonal	Dako	1:200
PAX5	24/Pax-5	BD Biosciences	1:100
CD3	LN10	Biocare Medical	1:200
Toluidine blue	Not applicable	Sigma	Not applicable
c-Kit	EP10	Biocare Medical	1:50
MUM1	BC5	Biocare Medical	1:50
CD20	Polyclonal	Invitrogen	1:50
CD79 α	HM47/A9	Biocare Medical	1:50
CD138	EP201	Bio SB	1:50